


Environmental Chambers and Rooms

Environmental chambers and rooms are specifically designed to provide uniform environmental conditions. They are ideal for incubation and for entomology experiments.




Knowledge Inside and Out

Expect to Grow

Accurate, easy to use controls drive successful research:

- Flexible shelving and lighting options are available
- Specifically suited for the constant temperature requirements of areas of discipline such as human tissue research (37C) and cold room storage (4C). Fully variable temperature ranges between -30C and 45C are also achievable
- Accurate and responsive temperature control
- Optional control of humidity, lighting, and CO₂ and other features



Environmental BioRoom™



Environmental Bigfoot™ Chamber

Expect to Know

Our BioChambers team will design and propose the best selection of chambers to maximize the growth area in your facility. Because we focus on solutions, you can expect growing success!



Our most popular series of products is found in research labs throughout the world. Bigfoot™ makes the best use of space. Bigfoot™ means more growth area....more experiments.....more results



BioBuild™ maximizes your growth area. Pack them together like dominos, side by side and back to back.



BioRoom™ enables scientists to conduct multiple experiments concurrently. Available in standard, specialized, and convertible rooms.



A state-of-the-art, independent controller providing accurate control of light, temperature, airflow, humidity, CO₂, watering and other options:

- Instantly shows you the status of all your experiments from your own computer
- Allows you to view and access any individual chamber with a click of the mouse



Available with every chamber and room is a unique proprietary monitoring system including the following features:

- A network of sensors which monitor both building services and the critical parameters of your experiment
- Multiple avenues of communication to ensure a prompt response to a cell phone, pager, PDA or similar device
- Remote, secure access—from anywhere

Engineering Excellence

BioChambers, formerly operating as Enconair, has been providing tailor-made, premium growth chambers and rooms to bio scientists around the world for over 30 years. At BioChambers, we're engineered to lead.

© Bio Chambers Incorporated, 2006. All rights reserved. Printed in Canada.

Knowledge Inside and Out

477 Jarvis Avenue
 Winnipeg, Manitoba
 Canada R2W 3A8
 Phone 204.589.8900
 Fax 204.582.1024
 www.biochambers.com